

Vials

- Septum Selection Guide
- Instrument Reference Chart
- Vial Size Quick Reference Chart
- Microsampling Vials & Inserts
- Crimp-Top Vials
- Screw-Thread Vials
- Certified Vials
- WISP Vials
- Headspace Vials
- Versa Vials
- VOA Sampling Vials
- Filter Vials, Sampling Valves for Vials
- Crimpers & Decappers



Septum Selection Guide

Materials	Compatibility	Incompatibility	Resealability	Max. Temp.
Red Rubber (synthetic)	acetone, alcohols, DMF, DMSO, ether	ACN, benzene, chloroform, heptane, hexane, pyridine, THF, toluene	very good	90 °C
PTFE/ Natural Rubber	PTFE: resistance until punctured Rubber: acetone, ACN, alcohols, diethylamine, DMF, DMSO, phenol	aromatics, carbon disulfide, chlorinated solvents, hydrocarbon solvents	very good	90 °C
PTFE/Silicone PTFE/Silicone/PTFE	PTFE: resistance until punctured Silicone: acetone, alcohols, DMF, DMSO, ether	ACN, benzene, chloroform, heptane, hexane, pyridine, THF, toluene	very good	205 °C
Polyethylene	Good resistance to solvents and weak acids or bases. Unreactive with most chemicals, but some solvents cause softening or swelling.	hydrocarbon solvents	one-time use	175 °C
Gray Chlorobutyl	acids or bases, water solutions, buffer solutions, oxygenated solvents, vegetable oils	aliphatic or aromatic hydrocarbons, halogenated solvents, mineral oils, strong acids	very good	100 °C

Abbreviations: ACN = acetonitrile, DMF = dimethylformamide, DMSO = dimethylsulfoxide, THF = tetrahydrofuran

NOTE: This chemical resistance chart is intended only as a guideline. It does not cover all compounds or all solvents. Tests were done at room temperature on pure, single solvents, and there is no data on solvent combinations. Always confirm the compatibility of your vial, closure, and chemical combination prior to sample preparations.

