


Organic Volatile Impurities (OVI) Analysis

G43 phase

Rxi®-624Sil MS Structure



Similar to: (6%-cyanopropyl(phenyl)-methyl)polysiloxane

similar phases

DB-624, VF-624ms, CP-Select 624 CB


Rxi®-624Sil MS Columns (fused silica)

(midpolarity Crossbond® phase)

- Low-bleed, high-thermal stability column—maximum temperatures up to 320 °C.
- Inert—excellent peak shape for a wide range of compounds.
- Selective—G43 phase highly selective for volatile organics and residual solvents, great choice for USP<467>.
- Manufactured for column-to-column reproducibility—well-suited for validated methods.

ID	df	temp. limits	20-Meter cat.#	30-Meter cat.#	60-Meter cat.#	75-Meter cat.#	105-Meter cat.#
0.18 mm	1.00 µm	-20 to 300/320 °C	13865				
0.25 mm	1.40 µm	-20 to 300/320 °C		13868	13869		
0.32 mm	1.80 µm	-20 to 300/320 °C		13870	13872		
0.53 mm	3.00 µm	-20 to 280/300 °C		13871	13873	13874	13875

Competitor Comparison: Class 2 - Mix A Residual Solvents for USP <467> Water-Soluble Articles



Improve system suitability pass rates with greater resolution on Rxi®-624Sil MS columns.

Column Rxi®-624Sil MS, 30 m, 0.32 mm ID, 1.80 µm (cat.# 13870)
Sample Residual solvents class 2 - mix A (cat.# 36271)
Diluent: Water
Injection Headspace-loop split (split ratio 5:1)
Liner: 1 mm split (cat.# 20972)
Headspace-Loop
Inj. Port Temp.: 140 °C
Instrument: Tekmar HT3
Inj. Time: 1 min
Transfer Line Temp.: 110 °C
Valve Oven Temp.: 110 °C
Sample Temp.: 80 °C
Sample Equil. Time: 60 min
Vial Pressure: 10 psi
Pressurize Time: 0.5 min
Pressure
Equilibration Time: 0.05 min
Loop Pressure: 5 psi
Loop Fill Time: 0.1 min
Oven
Oven Temp.: 40 °C (hold 20 min) to 240 °C at 10 °C/min (hold 20 min)
Carrier Gas He, constant flow
Linear Velocity: 35 cm/sec
Dead Time: 1.45 min @ 40 °C
Detector FID @ 250 °C
Data Rate: 5 Hz
Instrument Agilent/HP6890 GC
Acknowledgement Teledyne Tekmar

Peaks	tR (min)	Conc. (µg/mL)
1. Methanol	2.281	25.00
2. Acetonitrile	4.009	3.42
3. Dichloromethane	4.313	5.00
4. <i>trans</i> -1,2-Dichloroethene	4.798	7.83
5. <i>cis</i> -1,2-Dichloroethene	7.028	7.83
6. Tetrahydrofuran	7.706	5.75
7. Cyclohexane	8.708	32.33
8. Methylcyclohexane	14.099	9.83
9. 1,4-Dioxane	15.054	3.17
10. Toluene	22.018	7.42
11. Chlorobenzene	26.570	3.00
12. Ethylbenzene	26.837	3.07
13. <i>m</i> -Xylene	27.147	10.85
14. <i>p</i> -Xylene	27.147	2.53
15. <i>o</i> -Xylene	27.927	1.63

Organic Volatile Impurities (OVI) Analysis

Stabilwax® Columns (fused silica)


(polar phase; Crossbond® polyethylene glycol)

- Most stable polyethylene glycol (PEG) column available.
- Rugged enough to withstand repeated water injections.
- Lowest-bleed PEG column on the market; long column lifetimes.
- Temperature range: 40 °C to 260 °C.
- Equivalent to USP G14, G15, G16, G20, and G39 phases.

Restek's polar-deactivated surface tightly binds the Carbowax® polymer and increases thermal stability, relative to competitive columns. Because of the increased stability produced by the bonding process, Stabilwax® columns exhibit long column lifetimes, even when programming repeatedly up to 260 °C. The bonding mechanism of the column also produces polar compound retention times that do not shift, as is often observed on other wax-type columns. In addition, this bonding mechanism produces a column that can be rejuvenated by solvent washing.

ID	df	temp. limits	15-Meter cat.#	30-Meter cat.#	60-Meter cat.#
0.25 mm	0.10 µm	40 to 250/260 °C	10605	10608	10611
	0.25 µm	40 to 250/260 °C	10620	10623	10626
	0.50 µm	40 to 250/260 °C	10635	10638	10641
0.32 mm	0.25 µm	40 to 250/260 °C	10621	10624	10627
	0.50 µm	40 to 250/260 °C	10636	10639	10642
	1.00 µm	40 to 240/250 °C	10651	10654	10657
0.53 mm	0.25 µm	40 to 250/260 °C	10622	10625	10628
	0.50 µm	40 to 250/260 °C	10637	10640	10643
	1.00 µm	40 to 240/250 °C	10652	10655	10658
	1.50 µm	40 to 230/240 °C	10666	10669	10672
2.00 µm	40 to 220/230 °C	10667	10670		

ID	df	temp. limits	10-Meter cat.#	20-Meter cat.#
0.15 mm	0.15 µm	40 to 250/260 °C	43830	43831
0.18 mm	0.18 µm	40 to 250 °C		40602

G16 phase**Stabilwax® Structure****similar phases**

HP-INNOWax, CP-Wax 52 CB, VF-WAX MS, ZB-WAXplus

ordering note


Get the protection without the connection!
For Stabilwax® columns with built-in Integra-Guard® guard columns, see **page 23**.

free literature**Custom Residual Solvent Mixes**

Download your free copy from

www.restek.com

lit. cat.#
PHTS1212

**Residual Solvents (Class 1) on Stabilwax® (G16)**

Column Stabilwax®, 30 m, 0.32 mm ID, 0.25 µm (cat.# 10624)
Sample Residual solvents - class 1 (cat.# 36279)
 in 20 mL headspace vial (cat.# 24685)
 Water
Diluent:
Injection
 Inj. Vol.: 1.0 µL headspace-syringe split (split ratio 5:1)
 Liner: 2 mm splitless (cat.# 20712)
 Inj. Temp.: 140 °C
Headspace-Syringe
 Instrument: Overbrook Scientific HT200H
 Syringe Temp.: 100 °C
 Sample Temp.: 80 °C
 Sample Equil. Time: 45 min
 Inj. Speed: 8
 Inj. Dwell: 5 sec
Oven
 Oven Temp.: 50 °C (hold 20 min) to 165 °C at 6 °C/min (hold 20 min)
Carrier Gas He, constant flow
 Flow Rate: 2.15 mL/min
 Linear Velocity: 35.2 cm/sec
Detector FID @ 250 °C


Peaks
 1. 1,1-Dichloroethene
 2. 1,1,1-Trichloroethane
 3. Carbon tetrachloride
 4. Benzene
 5. 1,2-Dichloroethane

System suitability criteria met

GC_PH00951

G27 phase

Rtx®-5 Structure



Similar to: (5%-phenyl)-methylpolysiloxane

similar phases

HP-5, DB-5, CP-Sil 8 CB, ZB-5

NOTE: DB-5MS is a silarylene-based polymer, similar to Rxi-5Sil MS.

USP
Pharmaceutical
Standards

See pages 595–596.



Organic Volatile Impurities (OVI) Analysis

Rtx®-5 (G27) Columns (fused silica)

(low-polarity phase; Crossbond® diphenyl dimethyl polysiloxane)

- General-purpose columns for drugs, solvent impurities, pesticides, hydrocarbons, PCB congeners (e.g., Aroclor mixes), essential oils, semivolatiles.
- Temperature range: -60 °C to 350 °C.
- Equivalent to USP G27 and G36 phases.


The diphenyl dimethyl polysiloxane stationary phase is the most popular GC stationary phase and is used in a wide variety of applications. All residual catalysts and low molecular weight fragments are removed from the Rtx®-5 polymer, providing a tight mono-modal distribution and extremely low bleed.

ID	df	temp. limits*	15-Meter cat.#	30-Meter cat.#	60-Meter cat.#	105-Meter cat.#
0.25 mm	0.10 µm	-60 to 330/350 °C	10205	10208	10211	
	0.25 µm	-60 to 330/350 °C	10220	10223	10226	10229
	0.50 µm	-60 to 330/350 °C	10235	10238	10241	10244
	1.00 µm	-60 to 325/340 °C	10250	10253	10256	10259
0.32 mm	0.10 µm	-60 to 330/350 °C	10206	10209		
	0.25 µm	-60 to 330/350 °C	10221	10224	10227	
	0.50 µm	-60 to 330/350 °C	10236	10239	10242	
	1.00 µm	-60 to 325/340 °C	10251	10254	10257	10260
	1.50 µm	-60 to 310/330 °C	10266	10269	10272	10275
	3.00 µm	-60 to 280/300 °C	10281	10284	10287	10290
0.53 mm	0.10 µm	-60 to 320/340 °C	10207	10210		
	0.25 µm	-60 to 320/340 °C	10222	10225	10228	
	0.50 µm	-60 to 320/330 °C	10237	10240	10243	
	1.00 µm	-60 to 320/330 °C	10252	10255	10258	
	1.50 µm	-60 to 310/330 °C	10267	10270	10273	
	3.00 µm	-60 to 270/290 °C	10282	10285	10288	
	5.00 µm	-60 to 270/290 °C	10277	10279	10283	

ID	df	temp. limits	10-Meter cat.#	20-Meter cat.#	40-Meter cat.#
0.18 mm	0.20 µm	-60 to 325/340 °C	40201	40202	40203
	0.40 µm	-60 to 315/330 °C	40210	40211	40212

*Maximum temperatures listed are for shorter length columns. Longer columns may have a different maximum temperature.

Organic Volatile Impurities on Rtx®-5 (Rtx®-G27)



- Peaks**
1. Ethylene oxide
 2. Methanol
 3. Ethanol
 4. Diethyl ether
 5. 1,1-Dichloroethene
 6. Acetone
 7. Isopropanol
 8. Acetonitrile
 9. Methylene chloride
 10. *n*-Hexane
 11. *n*-Propanol
 12. Methyl ethyl ketone
 13. Ethyl acetate
 14. Tetrahydrofuran
 15. Chloroform
 16. 1,1,1-Trichloroethane
 17. Carbon tetrachloride
 18. Benzene
 19. 1,2-Dichloroethane
 20. Heptane
 21. Trichloroethylene
 22. *n*-Butanol
 23. 1,4-Dioxane
 24. Pyridine
 25. Toluene

Column Rtx®-5 w/5m Integra-Guard® Column (Rtx®-G27), 30 m, 0.53 mm ID, 5.00 µm (cat.# 10279-126)
Sample Headspace injection of common solvents for pharmaceutical processing. Prepared to equal about 500 ppm in the bulk pharmaceutical. Samples shaken and heated at 90 °C for 15 minutes, 1 mL headspace injection.

Injection
 Inj. Vol.: 1,000 µL headspace-syringe split (split ratio 2:1)
 Inj. Temp.: 220 °C

Oven
 Oven Temp.: 35 °C (hold 10 min) to 100 °C at 5 °C/min to 240 °C at 25 °C/min (hold 5 min)

Carrier Gas
 He, constant pressure
 Linear Velocity: 35 cm/sec @ 35 °C

Detector
 FID @ 240 °C

Notes
 FID sensitivity: 1.05 x 10⁻¹¹ AFS

Organic Volatile Impurities (OVI) Analysis

Rtx®-G27 Column (fused silica with 5-meter Integra-Guard® guard column)
(Crossbond® diphenyl dimethyl polysiloxane)

- Application-specific columns for residual solvents in pharmaceutical products.
- Analytical column with Integra-Guard® guard column eliminates connecting problems and leaks.
- Rtx®-G27 stable to 290 °C.

Some methods require the use of a guard column. Our Integra-Guard® integrated guard column system makes it easy to comply.

ID	df	temp. limits	30-Meter with 5-Meter, 0.53mm ID Integra-Guard Guard Column cat.#
0.53 mm	5.00 µm	-60 to 270/290 °C	10279-126

Rtx®-G43 Column (fused silica with 5-meter Integra-Guard® guard column)
(Crossbond® cyanopropylmethyl phenylmethyl polysiloxane)

- Application-specific columns for residual solvents in pharmaceutical products. Meet all requirements of USP <467>.
- Analytical column with Integra-Guard® guard column eliminates connecting problems and leaks.
- Rtx®-G43 stable to 240 °C.

Some USP <467> methods require the use of a guard column. Our Integra-Guard® integrated guard column system makes it easy to comply.

ID	df	temp. limits	30-Meter with 5-Meter, 0.53mm ID Integra-Guard Guard Column cat.#
0.53 mm	3.00 µm	-20 to 240 °C	16085-126

free literature

A Technical Guide for
Static Headspace
Analysis
Using GC

lit. cat.#
59895B



Custom Residual
Solvents
Mixes

lit. cat.#
PHTS1212




Download your free copies from

www.restek.com



USP <467> Residual Solvents on Rtx®-1301 (G43) by Static Headspace



Column

Rtx®-1301 w/5 m Integra-Guard®, 30 m, 0.53 mm ID, 3.00 µm
(cat.# 16085-126)

Sample

USP <467> calibration mixture #5 (cat.# 36007)

Diluent:

DMSO

Conc.:

To each 22 mL headspace vial 5ml water, ~ 1.0 g of sodium sulfate and 100 µL of stock standard were added. headspace-loop split (split ratio 2:1)

Injection

Headspace-Loop

180 °C

Inj. Port Temp.:

150 °C

Instrument:

Teledyne Tekmar HT3

Inj. Time:

1.0 min

Transfer Line Temp.:

150 °C

Valve Oven Temp.:

150 °C

Standby flow rate:

10 mL/min

Sample Temp.:

80 °C

Platen temp equil. time:

2.0 min

Sample Equil. Time:

15.0 min

Mixer time:

2.0 min

Mixing level:

5

Mixer stabilize time:

0.5 min

Vial Pressure:

15 psi

Pressurize Time:

2.0 min

Pressure Equilibration Time:

0.5 min

Loop Pressure:

5 psi

Loop Fill Time:

2.0 min

Loop fill equil. time:

0.5 min

Oven

40 °C (hold 20 min) to 240 °C at 25 °C/min (hold 10 min)

Carrier Gas

He, constant flow

Flow Rate:

5 mL/min

Detector

FID @ 250 °C

Make-up Gas Flow Rate:

45 mL/min

Notes

FID conditions:
hydrogen flow: 40 mL/min
air flow: 450 mL/min