

## Pinnacle® DB Columns: 1.9, 3, or 5 µm particles; 140 Å pore size

Restek® Pinnacle® DB columns are built for optimal UHPLC performance.

Pinnacle® DB columns are 100% manufactured by Restek in our Bellefonte, Pennsylvania, facility. Because performance begins with the support, our Pinnacle® DB UHPLC columns start with base-deactivated spherical silica that is optimized for UHPLC stability. From there, we bond them with a wide variety of phases to give chromatographers a stable and selective column. Get the most out of your UHPLC system. Combine selectivity and efficiency by using Restek® Pinnacle® DB UHPLC columns.




### Pinnacle® DB Aqueous C18 Columns (USP L1)

#### Chromatographic Properties

The Restek® Aqueous C18 is a rugged, reversed-phase column with a well-balanced retention profile. It can effectively retain more types of solutes than a conventional C18 and is ideal for multicomponent LC-MS analyses. The general-purpose Aqueous C18 boasts high reproducibility and compatibility with many mobile phase conditions—even 100% aqueous. And when used with a gradient, it eliminates the all-too-common issue of multiple compounds eluting near the column void time.

| Length                | 2.1 mm ID<br>cat.# | 3.0 mm ID<br>cat.# | 4.6 mm ID<br>cat.# |
|-----------------------|--------------------|--------------------|--------------------|
| <b>1.9 µm Columns</b> |                    |                    |                    |
| 30 mm                 | 9418232            | —                  | —                  |
| 50 mm                 | 9418252            | —                  | —                  |
| 100 mm                | 9418212            | —                  | —                  |
| <b>3 µm Columns</b>   |                    |                    |                    |
| 30 mm                 | 9418332            | 941833E            | 9418335            |
| 50 mm                 | 9418352            | 941835E            | 9418355            |
| 100 mm                | 9418312            | 941831E            | 9418315            |
| 150 mm                | 9418362            | 941836E            | 9418365            |
| <b>5 µm Columns</b>   |                    |                    |                    |
| 30 mm                 | 9418532            | 941853E            | 9418535            |
| 50 mm                 | 9418552            | 941855E            | 9418555            |
| 100 mm                | 9418512            | 941851E            | 9418515            |
| 150 mm                | 9418562            | 941856E            | 9418565            |
| 200 mm                | 9418522            | 941852E            | 9418525            |
| 250 mm                | 9418572            | 941857E            | 9418575            |




#### Column Characteristics:


|                    |  |
|--------------------|--|
| particle size:     | 1.9 µm, 3 µm, or 5 µm, spherical                       |
| pore size:         | 140 Å  |
| carbon load:       | 6%   |
| pH range:          | 2.5 to 8   |
| temperature limit: | 80 °C  |
| USP phase code:    | L1   |
| phase category:    | modified C18   |
| ligand type:       | proprietary polar modified and functionally bonded C18 |

#### Aqueous C18

#### USLC® Column Interaction Profile (See page 161 for more information.)



### Pesticides on Pinnacle® DB Aqueous C18 (LC-MS/MS, ESI+)



**Column** Pinnacle® DB Aqueous C18 (cat.# 9418252)  
**Dimensions:** 50 mm x 2.1 mm ID  
**Particle Size:** 1.9 µm  
**Pore Size:** 140 Å  
**Temp.:** 35 °C  
**Sample** multicomponent pesticide standard  
**Diluent:** water  
**Conc.:** 33.3 ppb each pesticide  
**Inj. Vol.:** 5 µL  
**Mobile Phase**  
**A:** 10 mM NH<sub>4</sub>OAc in water  
**B:** 10 mM NH<sub>4</sub>OAc in methanol

| Time (min) | %B |
|------------|----|
| 0          | 10 |
| 1          | 10 |
| 8          | 90 |
| 10         | 90 |
| 11         | 10 |

**Flow:** 0.60 mL/min  
**Max Pressure:** ~517 bar  
**Detector** Applied Biosystems/MDS Sciex LC-MS/MS  
**Model #:** 4000 QTRAP® LC-MS/MS system  
**Ion Source:** TurbolonSpray®  
**Ion Spray Voltage:** 5 kV  
**Gas 1:** 40 psi (275.8 kPa)  
**Gas 2:** 60 psi (413.7 kPa)  
**Source Temp.:** 500 °C  
**Instrument** Shimadzu UFLCxx




**Pinnacle® DB IBD UHPLC Columns (USP L68)**

**Chromatographic Properties**

The Restek® IBD is a polar-embedded column that acts as a strong hydrogen bonder and may be the most versatile column available today. With a unique polar group, this column is very retentive and selective for acids. It also provides symmetrical peak shape for strong bases. Restek's IBD is compatible with 100% aqueous mobile phases and can be used under reversed-phase or HILIC conditions to retain very polar, ionic compounds in highly organic mobile phases.

**Column Characteristics:**


|                    |   |
|--------------------|---|
| particle size:     | 1.9 µm, spherical                           |
| pore size:         | 140 Å                                       |
| end-cap:           | no  |
| pH range:          | 2.5 to 8                                    |
| temperature limit: | 80 °C                                       |
| USP phase code:    | L68   |
| phase category:    | polar-embedded alkyl                        |
| ligand type:       | proprietary polar functional embedded alkyl |



IBD

| Length                | 2.1 mm ID cat.# |
|-----------------------|-----------------|
| <b>1.9 µm Columns</b> |                 |
| 30 mm                 | 9425232         |
| 50 mm                 | 9425252         |
| 100 mm                | 9425212         |

**USLC® Column Interaction Profile**  
(See page 161 for more information.)



## Protect your column and your UHPLC performance with UltraShield and UltraLine UHPLC filters

A cost-effective way to extend the lifetime of any UHPLC column without sacrificing your UHPLC performance on any LC system.



See **page 188**.

[www.restek.com/LCguard](http://www.restek.com/LCguard)




## Pinnacle® DB Biphenyl Columns (USP L11)

### Chromatographic Properties

Since 2005, the Restek® Biphenyl has offered a greater degree of dispersion than conventional phenyls and a greater degree of polarizability than phenyl hexyls, creating higher selectivity and a greater range of usability. Because of these heightened interactions, this column shows substantial increases in retention—especially for dipolar, unsaturated, or conjugated solutes—and enhanced orthogonal selectivity when using methanol mobile phases. It is ideal for increasing sensitivity and selectivity in LC-MS analyses.

| Length                | 2.1 mm ID<br>cat.# | 3.0 mm ID<br>cat.# | 4.6 mm ID<br>cat.# |
|-----------------------|--------------------|--------------------|--------------------|
| <b>1.9 µm Columns</b> |                    |                    |                    |
| 30 mm                 | 9409232            | —                  | —                  |
| 50 mm                 | 9409252            | —                  | —                  |
| 100 mm                | 9409212            | —                  | —                  |
| <b>3 µm Columns</b>   |                    |                    |                    |
| 30 mm                 | 9409332            | 940933E            | 9409335            |
| 50 mm                 | 9409352            | 940935E            | 9409355            |
| 100 mm                | 9409312            | 940931E            | 9409315            |
| 150 mm                | 9409362            | 940936E            | 9409365            |
| <b>5 µm Columns</b>   |                    |                    |                    |
| 30 mm                 | 9409532            | 940953E            | 9409535            |
| 50 mm                 | 9409552            | 940955E            | 9409555            |
| 100 mm                | 9409512            | 940951E            | 9409515            |
| 150 mm                | 9409562            | 940956E            | 9409565            |
| 200 mm                | 9409522            | 940952E            | 9409525            |
| 250 mm                | 9409572            | 940957E            | 9409575            |




### Column Characteristics:


|                    |                                  |
|--------------------|----------------------------------|
| particle size:     | 1.9 µm, 3 µm, or 5 µm, spherical |
| pore size:         | 140 Å                            |
| carbon load:       | 8%                               |
| end-cap:           | yes                              |
| pH range:          | 2.5 to 8                         |
| temperature limit: | 80 °C                            |
| USP phase:         | L11                              |
| phase category:    | phenyl                           |
| ligand type:       | unique Biphenyl                  |

### Biphenyl

**USLC® Column Interaction Profile**  
(See page 161 for more information.)



## NSAIDs on Pinnacle® DB Biphenyl



- Peaks**
1. Uracil (void marker)
  2. Tenoxicam
  3. Sulfinyprazone

**Column** Pinnacle® DB Biphenyl (cat.# 9409565)  
**Dimensions:** 150 mm x 4.6 mm ID  
**Particle Size:** 5 µm  
**Pore Size:** 140 Å  
**Temp.:** 30 °C

**Sample**  
**Diluent:** 0.1% formic acid in water:methanol (40:60)  
**Conc.:** 100 µg/mL each component (see peak list)  
**Inj. Vol.:** 10 µL

**Mobile Phase**  
**A:** 0.1% formic acid in water  
**B:** methanol


| Time (min) | %B |
|------------|----|
| 0.00       | 60 |
| 2.0        | 60 |
| 8.0        | 90 |
| 20.0       | 90 |
| 20.1       | 60 |

**Flow:** 1.0 mL/min  
**Detector** UV/Vis @ 254 nm  
**Instrument** Shimadzu Prominence




**Column Characteristics:**

|                    |                                  |
|--------------------|----------------------------------|
| particle size:     | 1.9 µm, 3 µm, or 5 µm, spherical |
| pore size:         | 140 Å                            |
| carbon load:       | 6%                               |
| end-cap:           | yes                              |
| pH range:          | 2.5 to 8                         |
| temperature limit: | 80 °C                            |
| USP phase code:    | L43                              |
| phase category:    | fluorophenyl propyl              |
| ligand type:       | pentafluorophenyl propyl         |




PFP Propyl

**USLC® Column Interaction Profile**  
(See page 161 for more information.)



**Column Characteristics:**

|                    |                                  |
|--------------------|----------------------------------|
| particle size:     | 1.9 µm, 3 µm, or 5 µm, spherical |
| pore size:         | 140 Å                            |
| carbon load:       | 11%                              |
| end-cap:           | yes                              |
| pH range:          | 2.5 to 8                         |
| temperature limit: | 80 °C                            |
| USP phase code:    | L1                               |
| phase category:    | C18, octadecylsilane             |
| ligand type:       | monomeric C18                    |



C18

**Pinnacle® DB PFP Propyl Columns (USP L43)**

**Chromatographic Properties**

The Restek® PFP Propyl is a great choice for the retention and selectivity of charged bases, electronegative compounds, and amine-containing compounds. Unlike a conventional cyano column, the Restek® PFP Propyl is much more amenable to LC-MS because it is more reliable and efficient with acidic mobile phases. This versatile column is also compatible with highly aqueous mobile phases and HILIC separations.

| Length                | 2.1 mm ID cat.# | 3.0 mm ID cat.# | 4.6 mm ID cat.# |
|-----------------------|-----------------|-----------------|-----------------|
| <b>1.9 µm Columns</b> |                 |                 |                 |
| 30 mm                 | 9419232         | —               | —               |
| 50 mm                 | 9419252         | —               | —               |
| 100 mm                | 9419212         | —               | —               |
| <b>3 µm Columns</b>   |                 |                 |                 |
| 30 mm                 | 9419332         | 941933E         | 9419335         |
| 50 mm                 | 9419352         | 941935E         | 9419355         |
| 100 mm                | 9419312         | 941931E         | 9419315         |
| 150 mm                | 9419362         | 941936E         | 9419365         |
| <b>5 µm Columns</b>   |                 |                 |                 |
| 30 mm                 | 9419532         | 941953E         | 9419535         |
| 50 mm                 | 9419552         | 941955E         | 9419555         |
| 100 mm                | 9419512         | 941951E         | 9419515         |
| 150 mm                | 9419562         | 941956E         | 9419565         |
| 200 mm                | 9419522         | 941952E         | 9419525         |
| 250 mm                | 9419572         | 941957E         | 9419575         |



**Pinnacle® DB C18 Columns (USP L1)**

**Chromatographic Properties**

The general-purpose Restek® C18 is a conventional monomeric octadecylsilane column suitable for analyses of a wide range of compounds from acidic through slightly basic.

| Length                | 2.1 mm ID cat.# | 3.0 mm ID cat.# | 4.6 mm ID cat.# |
|-----------------------|-----------------|-----------------|-----------------|
| <b>1.9 µm Columns</b> |                 |                 |                 |
| 30 mm                 | 9414232         | —               | —               |
| 50 mm                 | 9414252         | —               | —               |
| 100 mm                | 9414212         | —               | —               |
| <b>3 µm Columns</b>   |                 |                 |                 |
| 30 mm                 | 9414332         | 941433E         | 9414335         |
| 50 mm                 | 9414352         | 941435E         | 9414355         |
| 100 mm                | 9414312         | 941431E         | 9414315         |
| <b>5 µm Columns</b>   |                 |                 |                 |
| 30 mm                 | 9414532         | 941453E         | 9414535         |
| 50 mm                 | 9414552         | 941455E         | 9414555         |
| 100 mm                | 9414512         | 941451E         | 9414515         |
| 150 mm                | 9414562         | 941456E         | 9414565         |
| 200 mm                | 9414522         | 941452E         | 9414525         |
| 250 mm                | 9414572         | 941457E         | 9414575         |

**also available**

**Trident Direct Guard Column System**

See page 189.



**Looking for an equivalent column?**

Restek has extensively studied column selectivity and can provide you with an accurate recommendation. Please contact Restek® Technical Support or your local Restek® representative.


## Pinnacle® DB C8 Columns (USP L7)



### Chromatographic Properties

Our C8 is a conventional monomeric octylsilane column offering a shorter alkyl chain to provide less hydrophobic retention and improved basic peak shape over a traditional C18 phase. Like our C18, this general-purpose Restek® C8 is suitable for a wide range of compounds from acidic through slightly basic.

| Length                | 2.1 mm ID<br>cat.# | 3.0 mm ID<br>cat.# | 4.6 mm ID<br>cat.# |
|-----------------------|--------------------|--------------------|--------------------|
| <b>1.9 µm Columns</b> |                    |                    |                    |
| 30 mm                 | 9413232            | —                  | —                  |
| 50 mm                 | 9413252            | —                  | —                  |
| 100 mm                | 9413212            | —                  | —                  |
| <b>3 µm Columns</b>   |                    |                    |                    |
| 30 mm                 | 9413332            | 941333E            | 9413335            |
| 50 mm                 | 9413352            | 941335E            | 9413355            |
| 100 mm                | 9413312            | 941331E            | 9413315            |
| <b>5 µm Columns</b>   |                    |                    |                    |
| 30 mm                 | 9413532            | 941353E            | 9413535            |
| 50 mm                 | 9413552            | 941355E            | 9413555            |
| 100 mm                | 9413512            | 941351E            | 9413515            |
| 150 mm                | 9413562            | 941356E            | 9413565            |
| 200 mm                | 9413522            | 941352E            | 9413525            |
| 250 mm                | 9413572            | 941357E            | 9413575            |



### Column Characteristics:

|                    |                                  |
|--------------------|----------------------------------|
| particle size:     | 1.9 µm, 3 µm, or 5 µm, spherical |
| pore size:         | 140 Å                            |
| carbon load:       | 6%                               |
| end-cap:           | yes                              |
| pH range:          | 2.5 to 8                         |
| temperature limit: | 80 °C                            |
| USP phase code:    | L7                               |
| phase category:    | C8, octylsilane                  |
| ligand type:       | monomeric C8                     |




## Pinnacle® DB Cyano Columns (USP L10)



### Chromatographic Properties

The Restek® Cyano is a traditional monomeric cyanopropylsilane that is recommended for assays where alternate selectivity, or confirmation, to a C18 or C8 column is desired. It can be used in normal-phase, reversed-phase (best with mobile phase pH between 5 and 7), and HILIC modes. It is an excellent choice for the analysis of protonated bases.

| Length                | 2.1 mm ID<br>cat.# | 3.0 mm ID<br>cat.# | 4.6 mm ID<br>cat.# |
|-----------------------|--------------------|--------------------|--------------------|
| <b>1.9 µm Columns</b> |                    |                    |                    |
| 30 mm                 | 9416232            | —                  | —                  |
| 50 mm                 | 9416252            | —                  | —                  |
| 100 mm                | 9416212            | —                  | —                  |
| <b>5 µm Columns</b>   |                    |                    |                    |
| 30 mm                 | 9416532            | 941653E            | 9416535            |
| 50 mm                 | 9416552            | 941655E            | 9416555            |
| 100 mm                | 9416512            | 941651E            | 9416515            |
| 150 mm                | 9416562            | 941656E            | 9416565            |
| 200 mm                | 9416522            | 941652E            | 9416525            |
| 250 mm                | 9416572            | 941657E            | 9416575            |



### Column Characteristics:

|                    |                           |
|--------------------|---------------------------|
| particle size:     | 1.9 µm or 5 µm, spherical |
| pore size:         | 140 Å                     |
| carbon load:       | 4%                        |
| end-cap:           | yes                       |
| pH range:          | 2.5 to 8                  |
| temperature limit: | 80 °C                     |
| USP phase code:    | L10                       |
| phase category:    | cyano                     |
| ligand type:       | cyanopropylsilane         |




Whether it's on our new Raptor™ SPP or the proven Pinnacle® DB and Ultra supports, Restek's LC Manufacturing group bonds our silica with stationary phases that offer maximum selectivity and reliability.



**Column Characteristics:**

|                    |                                  |
|--------------------|----------------------------------|
| particle size:     | 1.9 µm, 3 µm, or 5 µm, spherical |
| pore size:         | 140 Å                            |
| end-cap:           | no                               |
| pH range:          | 2.5 to 8                         |
| temperature limit: | 80 °C                            |
| USP phase code:    | L3                               |
| phase category:    | bare silica                      |
| ligand type:       | none                             |



Silica

**Pinnacle® DB Silica Columns (USP L3)**

**Chromatographic Properties**

Base-deactivated spherical silica is useful for normal-phase or HILIC separations.

| Length                | 2.1 mm ID cat.# | 3.0 mm ID cat.# | 4.6 mm ID cat.# |
|-----------------------|-----------------|-----------------|-----------------|
| <b>1.9 µm Columns</b> |                 |                 |                 |
| 30 mm                 | 9410232         | —               | —               |
| 50 mm                 | 9410252         | —               | —               |
| 100 mm                | 9410212         | —               | —               |
| <b>3 µm Columns</b>   |                 |                 |                 |
| 30 mm                 | 9410332         | 941033E         | 9410335         |
| 50 mm                 | 9410352         | 941035E         | 9410355         |
| 100 mm                | 9410312         | 941031E         | 9410315         |
| 150 mm                | 9410362         | 941036E         | 9410365         |
| <b>5 µm Columns</b>   |                 |                 |                 |
| 30 mm                 | 9410532         | 941053E         | 9410535         |
| 50 mm                 | 9410552         | 941055E         | 9410555         |
| 100 mm                | 9410512         | 941051E         | 9410515         |
| 150 mm                | 9410562         | 941056E         | 9410565         |
| 200 mm                | 9410522         | 941052E         | 9410525         |
| 250 mm                | 9410572         | 941057E         | 9410575         |



**Pinnacle® DB PAH UHPLC Columns**

**Chromatographic Properties**


Specifically designed to resolve complex mixtures of polycyclic aromatic hydrocarbons (PAHs). Get complete resolution of all 16 EPA 610 PAHs, plus two other routinely analyzed PAH compounds, in less than five minutes to greatly reduce run times and increase sample throughput.

**Column Characteristics:**

|                    |                   |
|--------------------|-------------------|
| particle size:     | 1.9 µm, spherical |
| pore size:         | 140 Å             |
| end-cap:           | no                |
| pH range:          | 2.5 to 8          |
| temperature limit: | 80 °C             |

| Length                | 2.1 mm ID cat.# |
|-----------------------|-----------------|
| <b>1.9 µm Columns</b> |                 |
| 30 mm                 | 9470232         |
| 50 mm                 | 9470252         |
| 100 mm                | 9470212         |

**PAH Mix on Pinnacle® DB PAH**



- Peaks**
1. Naphthalene
  2. Acenaphthylene
  3. 1-Methylnaphthalene
  4. 2-Methylnaphthalene
  5. Acenaphthene
  6. Fluorene
  7. Phenanthrene
  8. Anthracene
  9. Fluoranthene
  10. Pyrene
  11. Benzo[a]anthracene
  12. Chrysene
  13. Benzo[b]fluoranthene
  14. Benzo[k]fluoranthene
  15. Benzo[a]pyrene
  16. Dibenzo[a,h]anthracene
  17. Benzo[ghi]perylene
  18. Indeno[1,2,3-cd]pyrene

|                       |   |
|-----------------------|---|
| <b>Column</b>         | Pinnacle® DB PAH (cat.# 9470252)          |
| <b>Dimensions:</b>    | 50 mm x 2.1 mm ID                         |
| <b>Particle Size:</b> | 1.9 µm                                    |
| <b>Pore Size:</b>     | 140 Å                                     |
| <b>Temp.:</b>         | 30 °C                                     |
| <b>Sample</b>         | EPA Method 8310 PAH Mixture (cat.# 31841) |
| <b>Diluent:</b>       | acetonitrile                              |
| <b>Conc.:</b>         | 10 µg/mL                                  |
| <b>Inj. Vol.:</b>     | 1 µL                                      |

|                      |                                 |           |           |
|----------------------|---------------------------------|-----------|-----------|
| <b>Mobile Phase</b>  |                                 |           |           |
| A:                   | water                           |           |           |
| B:                   | acetonitrile                    |           |           |
| <b>Time (min)</b>    | <b>Flow (mL/min)</b>            | <b>%A</b> | <b>%B</b> |
| 0                    | 0.8                             | 60        | 40        |
| 2                    | 0.8                             | 40        | 60        |
| 4                    | 0.8                             | 0         | 100       |
| 4.5                  | 0.8                             | 0         | 100       |
| 4.51                 | 0.8                             | 60        | 40        |
| 5                    | 0.8                             | 60        | 40        |
| <b>Max Pressure:</b> | 724 bar                         |           |           |
| <b>Detector</b>      | Photo diode array @ 254, 4.8 nm |           |           |
| <b>Instrument</b>    | Waters                          |           |           |