

Sodium Sulfate (Bulk Adsorbent)

- Ideal for removing water from sample extracts.
- Activate by heating to 400 °C for four hours before use.
- Packaged in recloseable 5 kg buckets.

Anhydrous sodium sulfate is the most common drying agent used to remove moisture from sample extracts. We package our 60 mesh material in recloseable buckets.

Description	qty.	cat.#
Sodium Sulfate	5 kg	26204

**Florisil® PR** (Bulk Adsorbent)

- Pesticide residue grade.
- Packaged in glass containers.

Florisil® PR is commonly used to remove polar interferences from pesticide residues. This bulk material is ideal for labs packing their own extraction cartridges for pesticide residue extractions.

Description	qty.	cat.#
Florisil PR, 60/100 mesh	500 g	26135

**Granulated Activated Copper** (Bulk Adsorbent)

- Convenient form for removing sulfur from environmental extracts.
- Acidified and activated—ready for use.

Activated copper effectively removes elemental sulfur from environmental extracts. Our acid-washed and activated material can be used right out of the package. The 30 mesh granular material eliminates the potential for fine copper particles in filtered extracts.

Description	qty.	cat.#
Granulated Activated Copper, 30 mesh	1 kg	26136

**Ottawa Sand** (Bulk Adsorbent)

- Sample medium for matrix spikes and laboratory control blanks.
- Packaged in convenient 5 kg buckets.

Ottawa sand is listed in several U.S. EPA methods as the specified medium for matrix spike and laboratory control blanks.

Description	qty.	cat.#
Ottawa Sand	5 kg	26137

**Diatomaceous Earth** (Bulk Adsorbent)

- Improves extraction efficiency.
- Adsorbs moisture from samples.

Diatomaceous earth is used as a filter aid to improve extraction efficiency of densely packed soils, such as clays. By mixing the sample with diatomaceous earth, recoveries can be improved and excess moisture can be absorbed. Packaged in a convenient 1 kg quantity.

Description	Similar to Dionex Part #	qty.	cat.#
Diatomaceous Earth, 30/40 mesh	062819	1 kg	26033

