



24956

Autosampler Syringes for Bruker/Varian 8000 Series GCs

- Hamilton and SGE syringes are designed and tested to meet critical autosampler performance needs.
- Needle point styles are designed to withstand multiple, rapid injections through a septum.

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	cat.#	qty.	Restek cat.#
Bruker/Varian 8100/8200 Syringes								
10 µL	RN	26s	1.99"	5	701 Varian	202880	ea.	24956

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	cat.#	qty.	Restek cat.#
Bruker/Varian 8035, 8100, & 8200 Syringes								
10 µL*	R	25	53 mm	S/Hole	10R-VA8X-II	002924	ea.	24852

Bruker/Varian CP 8400, CP8410 Syringes

10 µL	F	23	50 mm	Cone	10F-VA8400-5/0.63C	002951	ea.	24923
10 µL	F	26	50 mm	Bevel	10F-VA8400-5/0.47	002950	ea.	21202
10 µL	R	26	50 mm	Cone	10R-VA8400-5/0-0.47C	002952	ea.	24930

*Gas-tight syringe.



24857

Replacement Needles for Bruker/Varian 8000 Series Autosampler Syringes

SGE

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	qty.	Restek cat.#
10 µL	23	53 mm	S/Hole	037779	2-pk.	24857
10 µL*	25	53 mm	S/Hole	037780	ea.	24858
10 µL**	25	53 mm	S/Hole	037777	ea.	24854
10 µL	25	53 mm	2	037776	2-pk.	24855

*Straight needle shaft, 0.20 mm ID needle, for use with viscous samples.

**Tapered needle shaft, 0.12 mm ID needle.