

Autosampler Syringes for PerkinElmer Autosystem and Clarus 500 GCs

- Hamilton and SGE syringes are designed and tested to meet critical autosampler performance needs.
- Needle point styles are designed to withstand multiple, rapid injections through a septum.

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	cat.#	qty.	Restek cat.#
5 µL	ASN	23	2.756"	3	75ASN/PE	88035	ea.	24953
5 µL	ASN	26	2.756"	3	75ASN/PE	88040	ea.	24952

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	cat.#	qty.	Restek cat.#
0.5 µL	R	23	70 mm	Cone	0.5BR-PE-0.63	000478	ea.	24811
0.5 µL	R	26	70 mm	Cone	0.5BR-PE-0.47	000475	ea.	24810
5 µL	F	23	70 mm	Cone	5F-PE-0.63	001954	ea.	24813
5 µL*	F	23	70 mm	Cone	5F-PE-GT-0.63	001957	ea.	24815
5 µL	F	26	70 mm	Cone	5F-PE-0.47	001953	ea.	24812
5 µL*	F	26	70 mm	Cone	5F-PE-GT-0.47	001955	ea.	24814
50 µL	F	23	70 mm	Cone	50F-PE-0.63	004670	ea.	24816

*Gas-tight.




24952

24812

Replacement Needle and Plunger Kit for PerkinElmer Autosampler Syringes

SGE

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	qty.	Restek cat.#
0.5 µL	23	70 mm	Cone	033765	kit	24818



24818

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) PTFE Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip