

Electronic Rechargeable Crimpers and Decappers

- Easy to use; comfortable grip.
- Hundreds of operations from one charge.
- Adjustable crimping force.
- Large battery allows for faster charge cycle.
- All kits include a universal plug set for the battery charger.
- For 11 mm and 20 mm caps.
- One-year warranty.

These electronic crimpers and decappers have a rechargeable 6.4-volt battery installed. The large battery allows for quick, powerful operation and a faster charge cycle, minimizing downtime. Charging the battery takes thirty minutes to two hours. The universal plug set includes U.S., European, UK, and Australian plug clips.

Description	Size	qty.	cat.#
Electronic Crimper	11 mm	kit	22358
Electronic Crimper	20 mm	kit	22359
Electronic Decapper	11 mm	kit	22360
Electronic Decapper	20 mm	kit	22361
Replacement Battery for Electronic Crimpers and Decappers		ea.	22362



Universal plug set for battery charger included in all kits.

Easy Grip Manual Crimpers and Decappers

Our Easy Grip crimpers and decappers feature tough, light, plastic bodies with curved handles that improve hand comfort during use. The bottom-pull handle allows for a steadier hold than top-push handles. An easily viewed adjustment knob is located at the head of the product. The clearly shown + and - symbols, along with directional arrows, simplify adjusting the intensity of the crimp/decap as needed. The knob is also designed to indicate when the crimp setting has been reached, when it fully contacts the body of the product.

Description	Size	qty.	cat.#
Crimper	11 mm	ea.	23396
Crimper	20 mm	ea.	23398
Decapper	11 mm	ea.	23397
Decapper	20 mm	ea.	23399



Safer than removing caps with pliers, knives, or screwdriver blades!

Restek Electronic Leak Detector

Protect your instrument
and your analytical column

www.restek.com/leakdetector