

Manufacturer Instrument/Model #	11 mm Crimp-Top	9 mm Short Cap Screw-Thread	Headspace	8 mm Screw-Thread	10 mm Screw-Thread	4 mL WISP 4.8 Style	1 mL WISP 96 Style
A.I.M. CPS-100, 200	x			x			
A.I. 42 Vial Tray	x						
A.I. Headspace (HS: 10, 20 mL)			x				
Agilent 1042, 1050, 1080, 1082, 1084, 1090, 1100, 5890, 6850, 6890, 7670A, 7671A, 7672, 7673A/B, 7683, 7890, 8042	x	x					
Agilent Headspace (HS: 6, 10, 20 mL, Flat)			x				
Alcott 738	x			x	x		
Alltech 570	x			x	x		
Alltech 580 (standard tray)	x						
Altex	x				x		
AMS 42	x						
Antek 736 Unisampler, 738	x						
ASC	x			x	x		
Beckman 501, 502, 507				x	x		
Bruker LC 51						x	
Bruker/Varian 8100/8200	x	x		x	x		
Bruker/Varian Marathon	x	x					
Bruker/Varian 8035, 8000 Series, 8055, 8085, 3800 GC				x			
Bruker/Varian 9100/9090/9095	x			x			
Bruker/Varian Genesis (HS: 10, 20 mL Rounded)			x				
Bruker/Varian/Rainin Dynamax AI-1A, AI-200, Dynamax AI-3	x						
Bruker/Varian/Rainin Dynamax AI-3AI-1A, AI-2W				x			
Carlo Erba/Fisons 42 vial tray AS 800	x						
Carlo Erba/Fisons AS 800		x					
Carlo Erba/Fisons AS 800, 8000				x			
Carlo Erba HS 500, HS 800 (HS: 6, 10, 20 mL, Flat)			x				
Carnegie CMA-250/200	x						
Chrompack CP 9000 GC Series	x						
Cueni Headspace (HS: 10, 20 mL)			x				
CTC CTC A105S	x						
CTC Headspace (HS: 20 mL, Flat)			x				
Dani ALS 39.80, 86.80	x	x		x			
Dani HSS 39.50, HSS 86.50 (HS: 10, 20 mL)			x				
Dani SPT 37.50 (HS: 20 mL)			x				
Delsi	x						
Dynatech 42 vial tray	x	x		x			
Dynatech LC2000		x					
Fisons 42 vial tray, AS800	x						
GBC	x						
Gerstel	x						
Gilson 231 XL, 232 XL, 233 XL, Aspec XL	x	x		x			
Gilson Asted XL	x			x			
Gynkotek Gina, others	x			x			
Hitachi AS-2000, AS-6000	x	x					
Hitachi S6551				x			
Hitachi L-7200				x		x	x
Hitachi L-7250				x		x	
IBM	x			x			
ICI other than LC 1600	x			x			
Infochroma	x			x			

Manufacturer Instrument/Model #	11 mm Crimp-Top	9 mm Short Cap Screw-Thread	Headspace	8 mm Screw-Thread	10 mm Screw-Thread	4 mL WISP 4.8 Style	1 mL WISP 96 Style
Jasco 851-AS, AS-950	x			x			
Jasco LC800 & 900 series					x		
Kipp	x			x			
Kontron 360	x			x			
Kontron 460				x			
LDC Marathon, Promis	x	x					
LDC other than 713	x			x			
L.E.A.P. Technologies CTC A105S	x						
L.E.A.P. Technologies CTC A200S	x	x		x	x		
L.E.A.P. Technologies HS 500 (HS: 10, 20 mL Rounded)			x				
Magnus Scientific	x			x			
PerkinElmer AI-1	x	x					
PerkinElmer Autosystem GC (HS: 10, 20 mL, Rounded)	x		x	x			
PerkinElmer F40, F45, HS 6, HS 40 (HS: 10, 20 mL, Rounded)			x				
PerkinElmer Integral 4000 (HS: 10, 20 mL, Rounded)	x		x				
PerkinElmer ISS 100, LC 600 42 vial tray	x	x					
PerkinElmer ISS 200 (HS: 6 mL)	x	x	x				
PerkinElmer 420/B, 4900	x			x			
Pharmacia LKB 2157-010	x			x			
Phillips 4247, 4710	x			x			
Phillips LC-XP	x			x		x	
Polymer-Labs GPC 110/210	x						
Precision Sampling GC111, GC 311, LC 241-60	x			x			
Pye LCXP	x						
S.G.E. M280D	x			x			
Sedere	x			x			
Shimadzu AOC 20i	x				x		
Shimadzu AOC-14/1400, AOC-17				x		x	
Shimadzu AOC-8B/9, SIL-6A, SIL-6B, -9A, -8A				x			
Shimadzu HSS-2B (HS: 27 mL)			x				
Shimadzu SIL-2AS						x	
Shimadzu SIL-10A, 10Ai, 10AxL	x	x		x	x		
Siemens AS 32, AS 200	x						
Spark Holland SPH 125	x	x					
Spark Holland Marathon, Promis	x	x		x			
Spark Holland Triathlon, Midas	x			x			
Spectra-Physics 8875, 8880	x	x		x			
Talbot	x			x			
Tekmar 7000, 7000/7050 (HS: 10, 20 mL, Rounded)			x				
Thermo Scientific TRACE GC 2000, AS2000	x	x					
TOA ICA5450				x			
Tosca 1	x						
TosoHaas TSK-6080, AS-8010, AS-8020	x			x		x	
Tracor 770, 771, 772	x						
TSP 8875, 8880, AS 100/1000, AS 300/3000	x			x			
Unicam 4710	x						
Unicam 4700LC				x			
Unicam 4247	x			x			
Unicam LC-XP	x			x		x	
United Technologies	x			x			
Waters Alliance 2690	x	x		x	x		
Waters 710, 717+						x	x
Waters 48-pos. M700						x	
Waters 96-pos. M700							x

