

Precleaned Volatile Organic Analyte (VOA) Sampling Vials

- Container, liner, and closure cleaned, assembled, and lot traceable.
- Open-top caps.
- PTFE faced 0.125" silicone septa.
- Cleaned to U.S. EPA Protocol B specifications.



24694



21798

26122

Description	Volume	Color	Material	Screw-Thread Size	qty.	cat.#
Precleaned VOA Vials	20 mL	Clear	Glass	24 mm/400	72-pk.	21798
Precleaned VOA Vials	20 mL	Amber	Glass	24 mm/400	72-pk.	21799
Precleaned VOA Vials	40 mL	Clear	Glass	24 mm/400	72-pk.	21796
Precleaned VOA Vials	40 mL	Amber	Glass	24 mm/400	72-pk.	21797
Collection Vials for ASE 200	60 mL	Clear	Glass	24 mm/400	72-pk.	26121
Collection Vials for ASE 200	60 mL	Amber	Glass	24 mm/400	72-pk.	26122
Replacement Septa, 24 mm x 0.125"			PTFE-lined silicone		100-pk.	24694



Restek® Certified Reference Materials

Crack open the best organic reference standards in the industry by choosing:

- ISO Quality
- Custom Flexibility
- Easy-to-Find Documentation

See **pages 464-465**.

www.restek.com/crm

