

**also available**

Limited Volume Inserts  
for 4.0 mL WISP 48 Screw-Thread  
Step Vials

See **page 356**.



**4.0 mL WISP 48 Crimp-Top 13 mm Vials, 15 x 45 mm, Snap Seal**



**4.0 mL WISP 48 Snap Seal Crimp-Top Vial** (vial only)

Description	Volume	Color	100-pk.	1,000-pk.
WISP 48 Snap Seal Vial w/White Graduated Marking Spot	4.0 mL	Clear	24658	24659



**Aluminum Seals w/Septa** for 4.0 mL WISP Crimp-Top Vials

Seal Color	Septa Material	100-pk.	1,000-pk.
Silver	PTFE/Butyl Rubber	21753	21754
Silver	PTFE/Silicone	21755	

**4.0 mL WISP 48 Screw-Thread Step Vials, 15 x 45 mm, 13/425**



**4.0 mL WISP 48 Screw-Thread Step Vials** (vials only)

Description	Volume	Material	100-pk.	1,000-pk.
WISP 48 Step Vial w/White Graduated Marking Spot	4.0 mL	Clear	24654	24655
WISP 48 Step Vial w/White Graduated Marking Spot	4.0 mL	Amber	24656	24657



Choose preslit caps to reduce the risk of needle bending, release vacuum from high-volume injections, and improve injection reproducibility when greater than 20% of vial volume is withdrawn.

**NEW!**




**4.0 mL WISP 48 Caps and Septa**

Description	Septa Material	100-pk.	1,000-pk.
WISP 48 Caps and Septa, Top Hat	PTFE/Silicone	21745	21746
WISP 48 Caps and Septa, Preassembled	Red PTFE/Silicone, 0.065"	21743	21744
Vial Caps & Septa, Black Polypropylene, Open-Hole	Red PTFE/Silicone, 0.060" w/Slit	26588	26589
WISP 48 Caps, Black Polypropylene, Open-Hole		21741	21742
Septa	Red PTFE/Silicone/PTFE, 0.040"	21751	21752
Septa	Red PTFE/Silicone, 0.060"	21749	21750

### 1.0 mL WISP 96 Style Shell Vials, 8 x 40 mm

Available in three styles.



also available

Limited Volume Inserts  
for 1.0 mL WISP 96 Style Shell Vials

See **page 356**.



### 1.0 mL WISP 96 Style, 8 x 40 mm Shell Vials

Description	Volume	Color	Material	200-pk.	1,000-pk.
WISP 96 Style Vial	1.0 mL	Clear	Glass	24662	24663
WISP 96 Style Vial	1.0 mL	Amber	Glass	24664	24665
WISP 96 Style Vial	1.0 mL	Clear	Polypropylene	24678	24679
WISP 96 Style Vial	700 µL Limited Volume	Clear	Polypropylene	24680	24681

8mm Starburst Plugs sold separately.




### 8 mm Starburst Plugs for Shell Vials, 1.0 mL WISP 96 Style, 8 x 40 mm

Description	Color	Material	200-pk.	1,000-pk.
Starburst Snap Plugs for Shell Vials	Clear	Polyethylene	21769	21770
Starburst Snap Plugs for Shell Vials, Conical	Clear	Polyethylene	21771	21772
Starburst Snap Plugs for Shell Vials, Silicone-Lined	Clear	Polyethylene	21773	21774



## Speed Up and Simplify GC Method Development with Restek's EZGC® Online Suite



Download these FREE web apps at [www.restek.com/ezgc](http://www.restek.com/ezgc)