



Toggle valve



Ball valve



Plug valve



23207



23209



23214



23211



23215

Shut-Off Gas Valves Swagelok®

Valve Type	1/8" Brass cat.#	1/4" Brass cat.#	1/8" Stainless Steel cat.#	1/4" Stainless Steel cat.#
Toggle	23142	23143	23198	23199
Ball	23144	23145	23200	23201
Plug	23146	23147	23202	23203

VCO® O-Ring Face Seal Fittings

(Stainless Steel & Treated) Swagelok®

- Unique design allows easy installation where space is limited.
- Assemblies can be used from high pressure to critical vacuum across a wide range of temperatures.
- Smooth finish on gland face ensures positive seal.
- Sealing is accomplished with a captive O-ring in the body component.

Swagelok® VCO® O-ring face seal fittings are designed for rapid assembly in pipe, tube, and welded systems.



Specifications:

Pressure Ratings	Up to 15,400 psig (1,061 bar)
Temperature Ratings	Up to 400 °F (204 °C)

Metering Gas Valves

Swagelok®



23206

Description	qty.	cat.#
1/8" Chrome-Plated Brass Metering Valve, straight	ea.	23148
1/4" Chrome-Plated Brass Metering Valve, straight	ea.	23149
1/8" Stainless Steel Metering Valve, straight	ea.	23204
1/4" Stainless Steel Metering Valve, straight	ea.	23205
Vernier Knob for Metering Valve	ea.	23206



23204



Toggle valve



Ball valve



Plug valve

Shut-Off Gas Valves Parker Balston®

Valve Type	1/8" Brass cat.#	1/4" Brass cat.#	1/8" Stainless Steel cat.#	1/4" Stainless Steel cat.#
Toggle	22188	22189	22190	22191
Ball	22192	22193	22194	22195
Plug	22196	22197	22198	22199

Description	Material	qty.	cat.#
1/4" VCO to 1/4" Tube Adaptor	Stainless Steel	2-pk.	23207
1/4" VCO to 1/4" Tube Adaptor	Siltek Treated	2-pk.	23209
1/4" VCO to 1/4" Tube Fitting	Siltek Treated	ea.	23211
1/4" VCO to 1/8" Tube Fitting	Stainless Steel	2-pk.	23215
1/4" VCO Body Blind	Stainless Steel	2-pk.	23214
Replacement O-Ring, 70 Durometer Fluorocarbon FKM, VCO 1/4", Size 010		10-pk.	23208

Metering Gas Valves

Parker Balston®



22209

Description	qty.	cat.#
1/8" Nickel-Plated Brass Metering, straight	ea.	22200
1/4" Nickel-Plated Brass Metering, straight	ea.	22201
1/8" Stainless Steel Metering, straight	ea.	22204
1/4" Stainless Steel Metering, straight	ea.	22205
Vernier Knob for Metering Valve	ea.	22209



22200