







Topaz Inlet Liners for Bruker/Varian GCs

COLUMN INSTALLS THIS END






Topaz Inlet Liners for Bruker/Varian GCs Equipped with 1177 Split/Splitless Inlets

	ID x OD x Length	Similar to Bruker/Varian Part #	qty.	cat.#
 4.0 mm ID Frit Liner	4.0 mm x 6.3 mm x 78.5 mm		5-pk.	23330
 4.0 mm ID Precision Liner w/Wool	4.0 mm x 6.3 mm x 78.5 mm		5-pk.	23328
 4.0 mm ID Single Taper Liner	4.0 mm x 6.5 mm x 78.5 mm	392611927	5-pk.	23331
 4.0 mm ID Single Taper Liner w/Wool	4.0 mm x 6.5 mm x 78.5 mm	392611936	5-pk.	23332

Topaz Inlet Liners for Bruker/Varian GCs Equipped with SPI Inlets

	ID x OD x Length	Similar to Bruker/Varian Part #	qty.	cat.#
 0.53 mm ID Straight Liner	0.53 mm x 4.6 mm x 54 mm	190010906	5-pk.	23460
 0.80 mm ID Straight Liner	0.80 mm x 4.6 mm x 54 mm	190010907	5-pk.	23461

Topaz Inlet Liners for Bruker/Varian GCs Equipped with 1078/1079 Inlets

	ID x OD x Length	Similar to Bruker/Varian Part #	qty.	cat.#
 3.4 mm ID Liner with Top Taper	3.4 mm x 5.0 mm x 54 mm	392611945	5-pk.	23329
 3.4 mm ID Frit Liner with Top Taper	3.4 mm x 5.0 mm x 54 mm	392611946	5-pk.	23462
 2.0 mm ID Liner with Top Taper	2.0 mm x 5.0 mm x 54 mm	392611947	5-pk.	23463
 0.75 mm ID Straight Liner	0.75 mm x 5.0 mm x 54 mm	392611948	5-pk.	23465
 3.4 mm ID Precision Top Taper Liner w/Wool	3.4 mm x 5.0 mm x 54 mm		5-pk.	23466

Patent pending



Topaz
INLET LINERS



Perfect Installation Every Time

There is never a question about which way to install a Topaz liner—simply orient it so that the column installs toward the “R” in the Restek logo.

