






Liners for APEX ProSep 800 & ProSep 800 Plus GCs

APEX Liners for ProSep 800 & ProSep 800 Plus GCs	Benefits/Uses:	ID* OD x Length	Similar to APEX Part #	ea. cat.#
 Mega IV (4.0 mm ID)	injections ≤125 µL	4.0 mm 6.0 mm x 243 mm	L-00410	21075
 Micro I (1.0 mm ID)	injections ≤5 µL	1.0 mm 6.0 mm x 243 mm	L-00110	21073
 MIDI II (2.0 mm ID)	injections ≤25 µL	2.0 mm 6.0 mm x 243 mm	L-00210	21074

*Nominal ID at syringe needle expulsion point.

Liners for ATAS Injectors

Liners for ATAS Injectors	Benefits/Uses:	ID* OD x Length	Similar to ATAS Part #	ea. cat.#	5-pk. cat.#
 ATAS Fritted Single Taper	dirty samples	3.0 mm 5.0 mm x 80 mm	A100126	22419	22420
 Liners for ATAS LINEX DMI System OPTIC 2 and 3 inlets		ID* OD x Length			5-pk. cat.#
ATAS LINEX DMI Liner		3.3 mm 5.0 mm x 81 mm			22353

*Nominal ID at syringe needle expulsion point.

COLUMN INSTALLS THIS END

Reference Standards Documentation Search

Search by cat.# or lot #!

- MSDSs
- Certificates
- Datapacks

www.restek.com/documents

