









Topaz Inlet Liners for DANI GCs

Topaz Inlet Liners for DANI GCs	ID x OD x Length	qty.	cat.#
 0.8 mm ID Straight Liner	0.8 mm x 6.0 mm x 70 mm	5-pk.	23293
0.8 mm ID Straight Liner (includes: ferrule set)	0.8 mm x 6.0 mm x 70 mm	5-pk.	23286
 2.0 mm ID Straight Liner	2.0 mm x 6.0 mm x 70 mm	5-pk.	23292
2.0 mm ID Straight Liner (includes: ferrule set)	2.0 mm x 6.0 mm x 70 mm	5-pk.	23285
 2.2 mm ID Straight Liner w/Wool	2.2 mm x 4.0 mm x 85.7 mm	5-pk.	23291
 2.2 mm ID Dimpled Liner w/Wool	2.2 mm x 4.0 mm x 55 mm	5-pk.	23296
2.2 mm ID Dimpled Liner w/Wool (includes: ferrule set)	2.2 mm x 4.0 mm x 55 mm	5-pk.	23289
 4.0 mm ID Straight Liner	4.0 mm x 6.0 mm x 70 mm	5-pk.	23294
4.0 mm ID Straight Liner (includes: ferrule set)	4.0 mm x 6.0 mm x 70 mm	5-pk.	23287
 4.0 mm ID Straight Liner w/Wool	4.0 mm x 6.0 mm x 70 mm	5-pk.	23295
4.0 mm ID Straight Liner w/Wool (includes: ferrule set)	4.0 mm x 6.0 mm x 70 mm	5-pk.	23288
 4.0 mm ID Single Taper Liner w/Wool	4.0 mm x 6.0 mm x 70 mm	5-pk.	23297
4.0 mm ID Single Taper Liner w/Wool (includes: ferrule set)	4.0 mm x 6.0 mm x 70 mm	5-pk.	23290

Topaz Inlet Liners for PerkinElmer GCs

Topaz Inlet Liners for PerkinElmer Auto SYS and Clarus 580/680 GCs Equipped with Split/Splitless Inlets: Split Injection	ID x OD x Length	Similar to PerkinElmer Part #	qty.	cat.#
 4.0 mm ID Straight Liner w/Wool	4.0 mm x 6.2 mm x 92.1 mm	N6502009	5-pk.	23449
 4.0 mm ID Precision Liner w/Wool	4.0 mm x 6.2 mm x 92.1 mm	N6121020	5-pk.	23450
Topaz Inlet Liners for PerkinElmer Auto SYS and Clarus 580/680 GCs Equipped with Split/Splitless Inlets: Splitless Injection	ID x OD x Length	Similar to PerkinElmer Part #	qty.	cat.#
 2.0 mm ID Straight Liner w/Wool	2.0 mm x 6.2 mm x 92.1 mm	N6121021	5-pk.	23451
Topaz Inlet Liners for PerkinElmer GCs Equipped with PSS Inlets	ID x OD x Length	Similar to PerkinElmer Part #	qty.	cat.#
 2.0 mm ID Straight Liner w/Wool	2.0 mm x 4.0 mm x 86.2 mm	N6121004	5-pk.	23317
Topaz Inlet Liners for PerkinElmer Auto SYS and Clarus 580/680 GCs Equipped with Split/Splitless Inlets: Direct Injection	ID x OD x Length	Similar to PerkinElmer Part #	qty.	cat.#
 4.0 mm ID Uniliner Liner w/Wool	4.0 mm x 6.2 mm x 92.1 mm	N6502016	5-pk.	23452
 4.0 mm ID Drilled Uniliner Liner (hole near top)	4.0 mm x 6.2 mm x 92.1 mm	N6121022	5-pk.	23453
 4.0 mm ID Drilled Uniliner Liner (hole near bottom)	4.0 mm x 6.2 mm x 92.1 mm	N6502013	5-pk.	23454
Topaz Inlet Liners for PerkinElmer GCs Equipped with PSS Inlets: Direct Injection	ID x OD x Length	qty.	cat.#	
 2.0 mm ID Drilled Uniliner (hole near top)	2.0 mm x 4.0 mm x 86.2 mm	5-pk.	23281	

Patent pending