


Replacement Cables


for Agilent GCs, Integrators, & Autosamplers

- Considerable savings over original manufacturer's equipment.
- Instructions and wiring diagrams included.
- Manufactured with only the highest-quality components.



Connect an Agilent 5890 GC to a non-Agilent integrator or standard strip chart recorder.

Description	Similar to Agilent Part #	length	qty.	cat.#
Replacement Cable	05890-60800	6 ft	ea.	20652



Connect an Agilent 5890 GC to remote-start another piece of equipment or to start the Agilent 5890 GC from that piece of equipment.

Description	Similar to Agilent Part #	length	qty.	cat.#
Replacement Cable	05890-61080	6 ft	ea.	20657



Injector Mounting Posts and Parking Post

for Agilent 7673 & 7683 Autosamplers

Performance equivalent to original manufacturer's parts.

Description	Similar to Agilent Part #	qty.	cat.#
Injector Mounting Post for Agilent 7673 Series Autosamplers for use with 6890 GCs	07673-21140	ea.	21237
Injector Mounting Post for Agilent 7683 Series Autosamplers for use with 6850/6890 GCs	G2613-20500	ea.	21172
Parking Post for Agilent 7673/7683 Series Autosamplers for use with 5890/6890/7890 GCs	05890-61525	ea.	22343

Turret Tray Assembly

for Agilent 7673 Autosamplers

Holds sample, waste, and solvent vials.



22855

Description	Similar to Agilent Part #	qty.	cat.#
Turret Tray Assembly for Agilent 7673 Autosamplers	07673-60605	ea.	22855

Autosampler Plunger Carrier Belt

for Agilent 7673A & 7673B Autosamplers

All parts meet or exceed original manufacturer's performance.



22695

Description	For Agilent Model #	qty.	cat.#
Autosampler Plunger Carrier Belt	7673A, 7673B	ea.	22695

Carriage Motor Belt

for Agilent 7673A, 7673B, and 7673C Autosamplers

Meets original manufacturer's performance.



22692

Description	Similar to Agilent Part #	qty.	cat.#
Carriage Motor Belt for Agilent 7673A, 7673B, and 7673C	1500-0676	ea.	22692

Z-Belt for Agilent 7673B & 7683 Autosamplers

Meets original manufacturer's performance.




22363

Description	Similar to Agilent Part #	qty.	cat.#
Z-Belt for Agilent 7673B & 7683 Autosamplers	1500-0803	ea.	22363

Injector Turret Motor for Agilent 7673A & 7673B Injectors

Meets original manufacturer's performance.



22337

Description	Similar to Agilent Part #	qty.	cat.#
Injector Turret Motor for Agilent 7673A & 7673B Injectors	07673-60810	ea.	22337