



Opti-Cap® GL-45 Bottle Tops

The most economical way to helium-sparge and deliver LC mobile phases. Opti-Cap® tops fit all standard GL-45 bottles and have two 1/8" holes and one 1/16" hole for tubing. All three openings are designed for threaded PEEK plugs.



25303

Description	qty.	cat.#
Opti-Cap (Cap and PEEK Plug)	ea.	25300
Opti-Cap Kit	kit	25301
Opti-Cap Kit with 1 L Bottle	kit	25302
Opti-Cap Kit with 2 L Bottle	kit	25303

Hub-Cap 4-Liter Bottle Tops

Hub-Cap bottle tops are a great way to neatly keep your mobile phase lines where they belong. Use them instead of plastic paraffin film, aluminum foil, or tape on your mobile phase reservoirs.



26541

Description	qty.	cat.#
Hub-Cap (assembly of the bottle cap and plug)	kit	26541
Hub-Cap Multi-Pack	3-pk.	26542

Hub-Cap Adaptors

Allow the use of the Opti-Cap® Top with 4-liter solvent bottles.

Description	qty.	cat.#
Hub-Cap Adaptor (allows use of the Opti-Cap w/4 L solvent bottles)	ea.	26538
Hub-Cap Adaptor Multi-Pack	3-pk.	26539
Hub-Cap Adaptor and Opti-Cap (for use w/GL-45 solvent bottles)	kit	26540

Related Items and Replacement Parts for Hub-Cap 4-Liter Bottle Tops

Description	qty.	cat.#
Mobile Phase Sparge Filter, 2 µm	ea.	25311
Mobile Phase Sparge Filter, 10 µm	ea.	25312
PTFE Tubing, 1/8" OD x 0.094" (2.4 mm) ID x 3 m	3 m	25307
PTFE Tubing, 1/8" OD x 0.063" (1.6 mm) ID x 3 m	3 m	25306
PEEK Plug, 1/4"-28 threads	3-pk.	25319

Opti-Cap® Adaptors

Allow the use of the Hub-Cap Top with GL-45 solvent bottles.



27197

26551

Description	qty.	cat.#
Opti-Cap Adaptor (allows use of Hub-Cap w/GL-45 solvent bottles)	ea.	27197
Opti-Cap Adaptor Multi-Pack	3-pk.	27198
Opti-Cap Adaptor w/Hub-Cap	kit	26551

Related Items and Replacement Parts for Opti-Cap® GL-45 Bottle Tops

Description	qty.	cat.#
Mobile Phase Sparge Filter, 2 µm	ea.	25311
Mobile Phase Sparge Filter, 10 µm	ea.	25312
PTFE Tubing, 1/8" OD x 0.094" (2.4 mm) ID x 3 m	3 m	25307
PTFE Tubing, 1/8" OD x 0.063" (1.6 mm) ID x 3 m	3 m	25306
PEEK Plug, 1/4"-28 threads	3-pk.	25319
1 L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25304
2 L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25305

Which Bottle Top Do I Need?

I have a 4-Liter bottle...

Option 1: Hub-Cap Top (cat.# 26541)

or

Option 2*: Opti-Cap® Top (cat.# 25300) plus Hub-Cap Adaptor (cat.# 26538)

*Also available as a kit (cat.# 26540)

Option 1



26541

Option 2*



25300

26538



*Also available as a kit (cat.# 26540)

I have a GL-45 bottle...

Option 1: Opti-Cap® Top (cat.# 25300)

or

Option 2*: Hub-Cap Top (cat.# 26541) plus Opti-Cap® Adaptor (cat.# 27197)

*Also available as a kit (cat.# 26551)

Option 1



25300

Option 2*



26541

27197



*Also available as a kit (cat.# 26551)