



25259



25270



25271



NEW!

25915

Preventive Maintenance (PM) Kits for Agilent HPLC Systems

Description	Model #	Similar to Agilent Part #	qty.	cat.#
Autosampler PM Kit Includes: rotor seal, needle seat, needle assembly, seat cap	1050	01078-68721	kit	25259
Autosampler PM Kit Includes: rotor seal, piston seals (2), needle assembly, needle seat, finger caps (3)	1100, 1200	G1313-68709	kit	25271
Pump PM Kit Includes: PTFE frits (2), outlet cap, gold disk seal, active inlet cartridge, piston seals (4), glass solvent filters (2)	1050, 1100, 1200	G1311-68710	kit	25270
Pump Maintenance Kit Includes: sapphire plunger (2), piston seals (2), outlet ball check valve, active inlet cartridge, PTFE frits (5)	1100, 1200	5065-4499	kit	25915



25258



25916, 25278



25272



25275




25349



26424

For even more options, visit www.restek.com/LCacc



25265

25058

Injector & Autosampler Supplies for Agilent HPLC Systems

Description	Model #	Similar to Agilent Part #	qty.	cat.#
Needle Seat	1050, 1090	79846-67101	ea.	25258
Needle Seat Assembly	1100, 1200, 1220	G1329-87017	ea.	25265
Needle Assembly	1100, 1200	G1313-87201	ea.	25278
Needle Assembly, 900 µL	1050, 1100, 1200	G1313-87202	ea.	25916
Rotor Seal, 2-Groove (not for use with 7125 injection valve)	1050	0101-0626	ea.	25272
Rotor Seal, 2-Groove	1100, 1200	0100-1853	ea.	25275
Rotor Seal (Rheodyne-Style), 3-Groove	1090	1535-4048	ea.	25349
Stator Face Assembly	1100, 1200	0100-1851	ea.	26424
Connecting Tube	1050, 1100	G1311-67304	ea.	25058
Injector Lubricant, PTFE, 10 g	1050, 1100	79841-65501	ea.	26520
Solvent Inlet Filter, SS	1050, 1090, 1100, 1200	01018-60025	ea.	26423
Isolation Seal for Injection Valve	1100, 1200	0100-1852	ea.	25917



26520



26423



25917


NEW!



26525




26530



26531




26526



26529



26532



26527



26533



26537



26528

Capillary Stainless Steel Tubing Assemblies for Agilent HPLC Systems

- Precut, micropolished tubing and preassembled fittings for quick, easy maintenance of your Agilent HPLC systems.
- Meet or exceed manufacturer's performance.

Description	Model #	Similar to	qty.	cat.#
		Agilent Part #		
Capillary SS Tubing With Fittings, 130 mm x 0.17 mm ID	1090	01090-87305	ea.	26525
Capillary SS Tubing With Fittings, 800 mm x 0.17 mm ID	1050	01078-87305	ea.	26526
Capillary SS Tubing With Fittings, 180 mm x 0.17 mm ID	1100	G1313-87305	ea.	26527
Capillary SS Tubing With Fittings, 700 mm x 0.25 mm ID	1050	01018-67305	ea.	26528
Capillary SS Tubing With Fittings, 700 mm x 0.25 mm ID	1050	01078-87306	ea.	26529
Seat Capillary, SS Tubing, 0.17 mm ID	1050	01078-87303	ea.	26530
Capillary SS Tubing, 105 mm x 0.17 mm ID	—	5021-1816	ea.	26531
Mixing Capillary Assembly	1100	G1312-67302	ea.	26532
Capillary SS Tubing, Valve to Metering Head	1100	G1313-87301	ea.	26533
Capillary SS Tubing, 150 mm x 0.17 mm ID	—	5021-1817	ea.	26534
Capillary SS Tubing, 280 mm x 0.17 mm ID	—	5021-1818	ea.	26535
Capillary SS Tubing, 400 mm x 0.17 mm ID	—	5021-1819	ea.	26536
1/16" Fitting, Front and Back Ferrules, Stainless Steel	—	5062-2418	10-pk.	26537



Pump Supplies for Agilent HPLC Systems

Description	Model #	Similar to Agilent Part #	qty.	cat.#
Piston Seals, PTFE w/Graphite, Black*	1050, 1100, 1200	5063-6589	2-pk.	22482
Piston Seals, PTFE w/Graphite, Black*	1050, 1100, 1200	5063-6589	10-pk.	22483
Piston Seals, Black, ASI	1050, 1100	5062-8516	10-pk.	25176
Piston Seals, Black	1090	5062-2494	4-pk.	25347
Pump Block Piston Seal	1050, 1100, 1200	0905-1420	2-pk.	25918
Seal Wash Kit, Binary Pump (4 seals, 4 gaskets)	1050, 1100, 1200	—	kit	25268
Seal Wash Kit (2 seals, 2 gaskets)	1050, 1100, 1200	—	kit	25269
Wash Seal	1050, 1100, 1200	0905-1175	ea.	25277
Active Inlet Cartridge	1050, 1100, 1200	5062-8562	ea.	26393
Seal, Gold Disk (outlet)	1050, 1100, 1200	5067-4728	ea.	25889
Outlet Cap	1050, 1100, 1200	5062-2485	4-pk.	25139
Outlet Cap & Gold Seal Assembly	1050, 1100, 1200	—	2-pk.	25890
Outlet Ball Valve, Binary Pump	1100, 1200	G1312-60012	ea.	25267
Outlet Ball Valve	1050, 1100, 1200	G1311-60012	ea.	25276
Sieves for Outlet Valve	1100, 1200	5063-6505	10-pk.	25266
Check Valve Cartridge Assembly	1090	79835-67101	ea.	25344
Frits, PTFE	1050, 1100, 1200	01018-22707	5-pk.	25466
Plunger Assembly	1050, 1100, 1200	5063-6586	ea.	25273
Sapphire Piston	1090	3980-0672	ea.	25345

*Graphite-filled PTFE; best for organic solvents.

For even more options, visit www.restek.com/LCacc



25261



25262



25399



25263



25888

Lamps for Agilent HPLC Systems

Description	Model #	Similar to Agilent Part #	qty.	cat.#
Detector Lamp, 1090 DA, 1050 VW/DA/MWD	1090, 1050	79883-60002	ea.	25260
Lamp, DAD G1315A, G1365A*	1100, 1200	2140-0590	ea.	25261
Lamp, VWD G1314A	1100, 1200	G1314-60100	ea.	25262
Lamp, DAD Long-Life Deuterium (2,000 hours)	1100, 1200 DAD	5181-1530	ea.	25399
8453 Deuterium Lamp	—	2140-0605	ea.	25263
Deuterium Lamp (2,000 hours)*	1100 DAD/MWD/1200	2140-0813	ea.	25888

*Compatible with 1260.

Other Supplies for Agilent HPLC Systems

Description	Model #	Similar to Agilent Part #	qty.	cat.#
Bottle Head Assembly	—	G1311-60003	ea.	25059