



Protect Your Column and Your UHPLC Performance With UltraShield and UltraLine UHPLC Filters

A cost-effective way to extend the lifetime of any UHPLC column without sacrificing your UHPLC performance on any LC system

Use with any UHPLC system



UltraShield UHPLC PreColumn Filter

- Cost-effective protection for UHPLC systems.
- Reliable way to extend column lifetime.
- Universal fit—connects easily to any brand column.
- Leak-tight to 15,000 psi (1,034 bar).
- 0.5 μm or 0.2 μm titanium frit in a stainless steel body with PEEK ferrule.

Specifications:

Inlet/Outlet:	Female/Male 10-32
Port Geometry:	Parker ($1/16$ CPI)
Material:	Titanium, stainless steel, PEEK ferrule
Filter:	0.5 μm or 0.2 μm stainless steel
Pressure Rating:	15,000 psig (1,034 bar)
Wrench Flat:	$5/16$ "

Description	Filter Porosity	qty.	cat.#
UltraShield UHPLC PreColumn Filter	0.5 μm frit	ea.	24995
		5-pk.	24996
		10-pk.	24997
UltraShield UHPLC PreColumn Filter	0.2 μm frit	ea.	25809
		5-pk.	25810
		10-pk.	25811

UltraLine UHPLC In-Line Filter

- In-line design installs easily with standard fittings.
- Cost-effective protection for UHPLC systems.
- Reliable way to extend column lifetime.
- Leak-tight to 15,000 psi (1,034 bar).
- Replaceable 0.5 μm stainless steel frit in stainless steel body.

Specifications:

Inlet/Outlet:	Female/Female 10-32
Port Geometry:	Parker ($1/16$ CPI)
Material:	Stainless steel housing
Filter:	0.5 μm stainless steel, 0.125" W x 0.062" T, 5 μL volume
Pressure Rating:	15,000 psig (1,034 bar)
Wrench Flat:	$3/8$ "

Description	qty.	cat.#
UltraLine UHPLC In-Line Filter (In-Line Assembly with Filter)	ea.	24993
UltraLine Replacement Filters	5-pk.	24994