

Atomic Absorption Calibration Standards

SCP SCIENCE has a wide selection of Standards; more than 31 elements available at 1000 µg/ml in 500 ml for Flame and Graphite Furnace Atomic Absorption

Satisfaction Guaranteed

"For best results, use SCP SCIENCE AA Single element Calibration Standards"

Calibration testing performed with NIST 3100 Series Standards

Certificate of Analysis with actual matrix, concentration and traceability to NIST


Long shelf life (up to 21 months unopened & 12 months opened)

Immediate availability

Free Standards with Lamps & GFAA Tubes*

AA Standard			Concentration 1000 µg/ml
		Matrix	Catalog Number
Aluminium	Al	HCl	140-002-135
Antimony	Sb	HNO ₃	140-001-515
Arsenic	As	HNO ₃	140-001-335
Barium	Ba	HNO ₃	140-001-565
Beryllium	Be	HNO ₃	140-001-045
Bismuth	Bi	HNO ₃	140-001-835
Boron	B	H ₂ O	140-000-055
Cadmium	Cd	HNO ₃	140-001-485
Calcium	Ca	HNO ₃	140-001-205
Chromium	Cr	HCl	140-002-245
Cobalt	Co	HNO ₃	140-001-275
Copper	Cu	HNO ₃	140-001-295
Gold	Au	HCl	140-002-795
Iron	Fe	HNO ₃	140-001-265
Lead	Pb	HNO ₃	140-001-825
Lithium	Li	HNO ₃	140-001-035
Magnesium	Mg	HNO ₃	140-001-125
Manganese	Mn	HNO ₃	140-001-255
Mercury	Hg	HNO ₃	140-001-805
Molybdenum	Mo	H ₂ O	140-000-425
Nickel	Ni	HNO ₃	140-001-285
Potassium	K	HNO ₃	140-001-195
Selenium	Se	HNO ₃	140-001-345
Silicon	Si	H ₂ O	140-000-145
Silver	Ag	HNO ₃	140-001-475
Sodium	Na	HNO ₃	140-001-115
Strontium	Sr	HNO ₃	140-001-385
Tin	Sn	HCl	140-002-505
Titanium	Ti	H ₂ O	140-000-225
Vanadium	V	HNO ₃	140-001-235
Zinc	Zn	HNO ₃	140-001-305

For a complete listing of **SCP SCIENCE** products, get your free Catalogue Today!!



www.scpscience.com

Sample Prep

ICP-AES & MS

AA

XRF

Calibration Standards

Certified Reference Materials

***Free Standards**



With every purchase of any Hollow Cathode Lamp or package of GFAA tubes, receive a 500 ml, 1000 µg/ml standard at no cost !!!

*Hazardous shipment charges are not included. The following element is excluded from all special offers: Au.

MK-MKGOO3-AAST-1.0-E

USA
Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499

Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67